

User Manual: Method for Evaluating Laboratory Accuracy of Infrared Thermographs

Tool Reference

RST Reference Number: RST24MD07.01

Date of Publication: 05/04/2023

Recommended Citation: U.S. Food and Drug Administration. (2023). *Infrared Thermographs-Methods for Evaluating Laboratory Performance* (RST24MD07.01). <https://cdrh-rst.fda.gov/infrared-thermographs-methods-evaluating-laboratory-performance>

For more information

[Catalog of Regulatory Science Tools to Help Assess New Medical Devices](#)

Disclaimer

About the Catalog of Regulatory Science Tools

The enclosed tool is part of the Catalog of Regulatory Science Tools, which provides a peer-reviewed resource for stakeholders to use where standards and qualified Medical Device Development Tools (MDDTs) do not yet exist. These tools do not replace FDA-recognized standards or MDDTs. This catalog collates a variety of regulatory science tools that the FDA's Center for Devices and Radiological Health's (CDRH) Office of Science and Engineering Labs (OSEL) developed. These tools use the most innovative science to support medical device development and patient access to safe and effective medical devices. If you are considering using a tool from this catalog in your marketing submissions, note that these tools have not been qualified as [Medical Device Development Tools](#) and the FDA has not evaluated the suitability of these tools within any specific context of use. You may [request feedback or meetings for medical device submissions](#) as part of the Q-Submission Program.

For more information about the Catalog of Regulatory Science Tools, email RST_CDRH@fda.hhs.gov.

Method for Evaluating Laboratory Accuracy of Infrared Thermographs

Acronym List

Acronym	Full name
ASTM	American Society for Testing and Materials
BB	blackbody
CS	calibration source (usually an extended area blackbody)
ETRS	external temperature reference source (usually a blackbody with known temperature and emissivity), often part of a ST
FOV	field of view
IEC	International Electrotechnical Commission (standards organization)
IR	Infrared
IRT	IR thermograph (also called IR camera), often works with an ETRS
ISO	International Organization for Standardization (standards organization)
MRTD	minimum resolvable temperature difference
NCIT	non-contact IR thermometers
ROI	region of interest
SD	standard deviation
ST	screening thermograph, often composed of an IRT and an ETRS
WTP	workable target plane (can be equal or less than the IRT's FOV)

Symbol List

Symbol	Definition
u_D	standard uncertainty of the measurement drift
u_S	standard uncertainty of the measurement stability
u_U	standard uncertainty of the measurement uniformity of the WTP
u_{ER}	standard uncertainty of the ETRS temperature
u_{MRTD}	standard uncertainty caused by the MRTD
u_{ST}	uncertainty of the ST $u_{ST}^2 = u_D^2 + u_S^2 + u_U^2 + u_{ER}^2 + u_{MRTD}^2$
u_{CS}	standard uncertainty of the CS used in calibration of the ST
u	combined standard uncertainty of the laboratory accuracy $u^2 = u_{CS}^2 + u_{ST}^2$
T_{ST}	temperature of the CS measured by the ST

T_{CS}	temperature of the CS
M_{8h}	mean value for 8 hours measurement
$M_{8h,min}$, $M_{8h,max}$	minimum and maximum of the mean value for 8 hours measurement
SD_{8h}	standard deviation for 8 hours of measurement

Setup, Device, and Environmental Requirements

- Minimum WTP resolution in thermograms: 320 × 240 pixels
- Minimum resolution for the subject's face in thermograms: 240 × 180 pixels
- Minimum ETRS resolution in thermograms: 20 × 20 pixels
- Minimum spatial resolution of ~1 mm /pixel (i.e., one pixel can image an area of 1×1mm² on the face)
- Ambient temperature: 18°C – 24°C
- Relative humidity: 10% – 75%
- Airflow from ventilation ducts should be deflected to minimize forced cooling or heating of the target.
- No source of IR radiation (for example, incandescent and halogen lighting) surrounds the experimental setup.

Comparison

[Our study](#) [2] evaluated and enhanced some requirements and test method in the IEC 80601-2-59:2017 standard [1]. The following table compares the differences between the standard methods and our proposed methods.

	Test procedure	<ul style="list-style-type: none"> • See clause 201.101.4. • The testing period should be the device's calibration interval or two weeks, whichever is longer. 	<ul style="list-style-type: none"> • Same as the standard. Detailed instruction is provided in Section 2.3.1. • The testing period should be the device's calibration interval or one month, whichever is longer.
	Requirements	Stability: <ul style="list-style-type: none"> • $3 \cdot SD_{8h} < 0.1^{\circ}\text{C}$ Drift: <ul style="list-style-type: none"> • $M_{8h,max} - M_{8h,min} < 0.1^{\circ}\text{C}$ The combined stability and drift: <ul style="list-style-type: none"> • Stability + Drift $< 0.2^{\circ}\text{C}$. 	Stability: <ul style="list-style-type: none"> • Standard uncertainty $u_s = SD_{8h} < 0.03^{\circ}\text{C}$, or • Expanded uncertainty $U_s = k \cdot u_s < 0.1^{\circ}\text{C}$, with coverage factor $k=3$ and confidence level $>99\%$. Drift: <ul style="list-style-type: none"> • Standard uncertainty $u_d = M_{8h,max} - M_{8h,min} < 0.1^{\circ}\text{C}$
	Notes	<ul style="list-style-type: none"> • Following the standard, we calculated the mean (M_{frame}) of the ROI on each captured frame, followed by the mean (M_{8h}) and SD (SD_{8h}) of all the M_{frame} values within one day (8 hours). $M_{8h,max}$ and $M_{8h,min}$ are the maximum and minimum M_{8h} values during the drift testing period. • Since the drift trend can be monotonical, drift is the maximum difference between M_{8h} values during the testing period (i.e., the device's calibration interval or two weeks, whichever is longer). • Calculation of M_{frame} will reduce random noise by a factor of \sqrt{N}, where N is the number of pixels in the ROI for mean calculation. The random noise is considered in the uniformity test. • Based on our recommendations, the requirement for combined stability and drift is redundant and can be removed. 	

	Test procedure	<ul style="list-style-type: none"> • See Clause 201.101.6. • The CS is repositioned at no less than 29 locations (four corners, center, and at least 24 random locations) throughout the WTP, and the CS temperature was measured with the IRT at each location. 	<ul style="list-style-type: none"> • See Section 2.3.2 for details. • The CS is repositioned in front of the IRT so that the CS aperture can cover the entire FOV of the IRT. • The camera should be focused on the working distance for uniformity evaluation to avoid the effect of focal length. • The SD of temperature values from all the pixels within the WTP is defined as the standard uncertainty of uniformity (u_U).
	Requirements	<ul style="list-style-type: none"> • The maximum difference between the 29 CS temperature values is less than 0.2°C. 	$u_U < 0.05^\circ\text{C}$
	Notes	<ul style="list-style-type: none"> • The standard does not explain whether the CS temperature should be read from one pixel or the average of several pixels. Since only 29 values are used, the test result has poor repeatability and large variation. 	
	Test procedure	<ul style="list-style-type: none"> • See Clause 201.101.5. Complies with the ASTM E1213-14 standard [4] • A standard four-bar target is mounted in front of the CS (a differential blackbody). • The target is imaged with the IRT and displayed on a video monitor where the image is viewed by an observer. • The temperature difference (ΔT) between the bars and their conjugates, initially zero, is increased incrementally only until the observer can distinguish the four bars. This critical temperature 	<ul style="list-style-type: none"> • See Section 2.3.3 for details. • Same setup with the ASTM E1213-14 method. • Instead of using observers, the contrast of the target bar images is digitally analyzed directly. The minimum visible contrast is defined as 10% and the ΔT value corresponding to the 10% contrast is defined as MRTD. • MRTD can be affected by the target/ambient temperature. Since human face temperature usually ranges from 34 °C to 36 °C, we recommend the bar target temperature for MRTD measuring is set within this range. • The MRTD values in the horizontal and vertical directions might be different and should be measured in both directions.

		<p>difference is the MRTD for this observer.</p> <ul style="list-style-type: none"> Repeat the procedure for multiple observers and report the median MRTD values for all the observer as the final MRTD result. The ASTM E1213-14 standard has no requirement for the spatial frequency and background temperature of the bar target, but requires these values be reported. Monitor quality and setting (e.g., brightness and contrast) and viewing environment might affect results. 	<ul style="list-style-type: none"> MRTD will increase with the bar target spatial frequency and should be evaluated with different bar targets of different spatial frequencies. A graph of MRTD versus spatial frequency is a convenient form of reporting the data. The highest spatial frequency should at least be 0.2 cycles/mrad. The MRTD uncertainty (u_{MRTD}) is defined as the difference of MRTD values between the highest and lowest target spatial frequencies. The final u_{MRTD} value is the average of u_{MRTD} values in both the vertical and horizontal directions.
	Requirements	MRTD is not larger than 0.1°C	MRTD at the highest spatial frequency (≥ 0.2 cycles/mrad) is not greater than 0.1°C.
	Notes	<ul style="list-style-type: none"> The ASTM E1213-14 method is subjective, and thus susceptible to reader variability, which is not ideal for standardization. By omitting observers, our proposed method significantly simplifies the MRTD measurement. The ASTM E1213-14 method does not specify MRTD value at which spatial frequency should satisfy the requirement. The method does not discuss the effect of target temperature and direction either. 	
	Requirements	<ul style="list-style-type: none"> The measured radiometric temperature laboratory accuracy shall satisfy $T_{ST} - T_{CS} + u \leq 0.5$ over the range of at least 34°C to 39°C at no less than five CS temperature points. 	
	Requirements	≥ 0.998 (Annex BB)	0.98 ± 0.01
	Note	It is difficult to find a commercial blackbody with emissivity ≥ 0.998 in the temperature range of 30°C to 40°C (clause 201.101.2.1). While some cavity blackbodies have emissivity around 0.99, they usually operate at temperature higher than	

		50°C. We were unable to identify a justification for the specified emissivity in the standard (references in the standard do not mention an emissivity of 0.998). We believe CS emissivity of 0.98 ± 0.01 should be sufficient if proper emissivity compensation is applied in the IRT algorithm.	
	Requirements	Less than 10% of the face (Annex AA)	A larger size ($15 \pm 20\%$) is also acceptable if it can be experimentally proven that the ETRS does not adversely affect the measurement.

Note: Contents in this appendix are based on references [1, 2].

References:

- [1] IEC 80601-2-59: [Medical electrical equipment - Part 2-59: Particular requirements for the basic safety and essential performance of screening thermographs for human febrile temperature screening](#), IEC, Geneva, Switzerland, 2017.
- [2] P. Ghassemi, T. J. Pfefer, J. P. Casamento, R. Simpson, and Q. Wang, "[Best practices for standardized performance testing of infrared thermographs intended for fever screening](#)," *PLoS ONE*, vol. 13, no. 9, p. e0203302, 2018.
- [3] Q. Wang, Y. Zhou, P. Ghassemi, D. McBride, J. P. Casamento, and T. J. Pfefer, "[Infrared Thermography for Measuring Elevated Body Temperature: Clinical Accuracy, Calibration, and Evaluation](#)," *Sensors*, vol. 22, p. 215, 2022.
- [4] ASTM E1213-14: [Standard Practice for Minimum Resolvable Temperature Difference for Thermal Imaging Systems](#), ASTM, West Conshohocken, PA, USA, 2014.